

have the Gorgeous, Full, Kissable Lips you always wanted... NATURALLY

Luscious Lips

www.lippump.com

HOW IT WORKS

The Luscious Lips device creates a vacuum that gently draws fluid into the lips. In the beginning, your lips are not accustomed to holding the additional fluid, so it is necessary to gradually build up the tissue by using the device gently and sparingly, especially in the critical first few days.

BEGINNERS: Condition your lips with the Luscious Lips device for the first 10-14 days.

Tilt your head back slightly and place the Luscious Lips device over your mouth, ensuring that the device is covering the entire mouth. Hold the device firmly in place and create an airtight seal by placing your thumb and forefinger on each side of the mouthpiece. To create an initial vacuum, suck the air out of the mouthpiece. You should feel the mouthpiece snug up to your face. (If you have a narrow face, it may be necessary to turn the device upside down to create an airtight seal).



Pull the outer cylinder slowly with the opposite hand and allow the device to slowly pull the lips gently into the cylinder. It is not necessary to pull the cylinder all the way out. Hold this position for TWO SECONDS and release the vacuum by returning the outer cylinder to its original position. Repeat this process four additional times for a total of five 2-second pulls each day for the first 10 to 14 days.

INTERMEDIATE LEVEL: Creating plump, full lips.

Following the initial conditioning period, increase the duration of the "pull" gradually each day, beginning with three 5-second pulls for two or three days and increasing to three 10-second pulls by the end of the first week following the initial conditioning period. You will see your lips plump before your very eyes! Three 10-second plumps will create a light plump. You may choose to re-plump two or three times a day.

ADVANCED LEVEL: Create the Luscious Lips you want.

As an intermediate user, you know how effective the Luscious Lips device can be, so use your own good judgment when creating the plumpness you desire in your lips. Most users seek to achieve a moderate plump with three or four 15-second plumps which lasts for a few hours. Then, they re-plump their lips two or three times a day to maintain the moderate plump. Do not exceed 120 seconds of use per day.

USE & CARE OF YOUR DEVICE:

Before you begin, lubricate your Luscious Lips device.



For smooth operation of your Luscious Lips device, we recommend that you occasionally lubricate the inside of the tube connected to the mouthpiece. It's easy to do. Simply pull the device apart and apply a small amount of facial moisturizing cream up INSIDE the tube connected to the mouthpiece. Reassemble the device and pump it a few times to create a smooth, easy pull. Now you're ready to begin conditioning your lips.

Cleaning Instructions:

Please do not immerse the Luscious Lips device in water. To clean, wipe with a tissue or soft cloth dampened, if necessary, with alcohol or a mild detergent.



get
Beautiful,
Smooth, Fuller
lips today!

call 800 592-3366

*"Luscious Lips
works for me and I
recommend it to anyone
who wants fuller,
younger looking lips."*

Carolyn Doherty, MD